



X-Pack Flat upgrade

Retrofit kit

Product Highlights

- For easy upgrade of analogue systems to digital
- Integrated with different X-ray Generators
- Bucky unit with removable grid
- Grid for optimal scatter reduction
- 5-Field Ionization Measuring chamber (measuring fields not visible)
- Detector with best in class image quality incl. image processing

Technical Data

5-Field Chamber

- Exposure time range 1ms – 6s
- Ionization voltage 360V – 550V
- Reproducibility 1%

X-Flat 17/17 S

- Pixel pitch 143um
- Time to Image 3,5 sec. after X-ray exposure
- Resolution 3,5 lp/mm

X-ray Generators

<u>Manufacturer</u>	<u>Types</u>
Sedecal, Spain	different models
Stadler, Switzerland	8X series
Control-X, Hungary	different models
CPI, Canada	Indico100, CMP (wip)
Shimadzu, Japan	different models (wip)
Röntgenwerk Bochum, Germany	different models (wip)
EMD, Canada	different models (wip)
PSM, Italy	different models (wip)

Bucky unit for 17" Detector

