



X-Tray

Bucky Tray for X-Ray Detectors

Product Highlights

- For easy upgrade of analogue systems to digital
- Grids can be changed or removed in few seconds
- Can be equipped with automated exposure chambers
- Suitable for all standard fix detectors
- Oscillating grid for better image quality and less dose
- Continuous ready mode for serial examination


Technical Data

Suitable Flat Panel Detectors

Vendor	Dimensions in mm		
	Width	Length	Height
Samsung	500	497	45
Trixel	500	495	50
Varian	469	469	36
Toshiba	512	495	43
DR Tech	482	482	35
Atlaim	500	500	40

X-Tray Parameters

- Dimension: 525 x 586 mm, Height 77 mm
- Weight: 8,2 kg



All rights are reserved. Reproduction in whole or parts is prohibited without the prior written consent of copyright holder.
X-Alliance reserves the right to make changes in specification and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Distributed by:
X-Alliance GmbH
Suhrenkamp 59-69
22335 Hamburg - Germany
Tel: +49 40 600 88 44 0
Fax: +49 40 600 88 44 44
www.x-alliance.com
sales@x-alliance.com