

Bucky WallStand	Table	Collimator/Tube	Generator, MEX
		XA ○ (6mm ²)	○ ⊕ MEX 106
		Ctrl Handle ○ (6mm ²)	○ ⊕ MEX 107
	Table PE1 ○	(8mm ²)	○ ⊕ MEX 108
	Bucky unit ○	(4mm ²)	○ ⊕ MEX 109
VP-Box ○		(8mm ²)	○ ⊕ MEX 110
Bucky unit ○ (6mm ²)			

Note: The ⊕ Symbol is indicated on the cables as "PE" (Protective Earth)

Bucky System

Earthing Plan

2012.1/4535-999-01241 Z6-1.2

© X-Alliance GmbH