



Sig. LED	Function	Color	Explanation
H1		green	+24V Supply
H2		green	Cass present Table
H3		green	Cass present Wallstand
H4		yellow	RQ Xray Table
H5		yellow	RQ Xray Wallstand
H6		yellow	Field R Table Chamber
H7		yellow	Field C Table Chamber
H8		yellow	Field L Table Chamber
H9		yellow	Field R Wall Chamber
H10		yellow	Field C Wall Chamber
H11		yellow	Field L Wall Chamber
H12		green	VCC +15V
H13		green	VEE -15V

M.P.	Function	Explanation
MP1	GND	electronic Ground
MP2	+15V	VCC chamber Table
MP3	0-12V Ramp	Chamber Table Dose Signal
MP4	+15V	VCC chamber Wall
MP5	0-12V Ramp	Chamber Wall Dose Signal
MP6	+400V	Chamber HV
MP7	Reset	Chamber Wall
MP8	Reset	Chamber Table Dose Signal
MP9	low active	Synch Table
MP10	low active	Synch Wallstand

W6	2-3	Service: Without Cass present Table
W7	2-3	Service: Without Cass present Wall

Bucky System

Location Plan

2012.1/4535-999-01011 Z3-1.0

© X-Alliance GmbH