



FLXIS

Dynamic X-Ray Image Intensifiers

Product Highlights

- Complete image chain and full range of innovative products covering all X-ray applications
- Customisation based on standard Philips components
- Decrease on-site engineering and integration costs
- Built in Philips medical application knowledge
- Extension to integrated image processing
- One stop shopping

Technical Features


- 9" / 12" / 16" II assembly
- Leading technology in image intensifiers
- Reliable, high resolution camera
- Integrated dose control
- Flexible interfacing & control

Technical Data / Configuration options

FLXIS Configuration Options	
9" / 12" II Assembly	
Camera	0,5k2 / 1k2 Fluoro camera
	0,5k2 / 1k2 Exposure camera
Subsystem control	
Image processing	Standard IP (LIH & ABC)
	IP Surgery
	IP Fluoro / IP Vascular
	IP Overlay / IP Plus

Options

- Dual monitor output (option to fluoroscopy package or surgery package)
- DVI video output / RS644 video output
- DVLP Gigabit Ethernet
- Direct Control Interface (DCM)
- X-ray grid (tbd specified)
- Mechanical fixation (front or lateral)
- Color RAL 9002 (other colors on request)



All rights are reserved. Reproduction in whole or parts is prohibited without the prior written consent of copyright holder.
X-Alliance reserves the right to make changes in specification and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Distributed by:
X-Alliance GmbH
Suhrenkamp 59-69
22335 Hamburg - Germany
Tel: +49 40 600 88 44 0
Fax: +49 40 600 88 44 44
www.x-alliance.com
sales@x-alliance.com