



## System Integration Support

### Competence Centre for X-Ray Imaging

- X-Alliance supports the integration of any components for x-ray diagnostics.
- Based on Philips know-how we provide highest quality imaging and application excellence.

### Co-Development

- X-Alliance has a number of exclusive distribution agreements for first-class OEM components and systems. We have own products like collimators and high voltage cable and the unique line of X-Packs, customized packages, pre-configured and adjusted.
- Our know-how in X-ray allows us to offer customization and co-development projects in the entire range of X-ray systems. Our team has many years of experience and supports the development of complete systems, subsystems, or upgrades - whatever you are just planning.

### Customization with Tailor-made Solutions

- We help you creating your system. Geometries from our partners with our packages and our development support enable you to develop analogue or digital solutions.
- Create your own customized system: components, colours and design according to your taste with our support.

## Time-to-market solutions

- Our proven packages enable you to provide solutions within short development times.
- We support your time-to-market, your system quality, and the unique look of your system: Every level of your development process: specification, design, engineering, documentation, labelling and certification.
- This gives you the freedom to develop new systems without having fixed headcounts in your organization.



## Examples

- **X-Pack** is a complete package of components for an X-ray system preconfigured and adjusted. Not only the interfacing of all components but also the image quality is fine-tuned to perfection for your application. X-Pack is available for Bucky systems, RF-systems, surgery systems as well as special applications (e.g. Lithotripsy systems) from 3.5 kW to 80 kW and from 9" Image Intensifiers to Flat Panel Detectors.
- **X-Pack Flat Upgrade** is a retrofit kit for easy upgrade of an analogue system to digital.
- **Bucky|Röntgen** is a readymade system solution, including design, documentation and certification.
- **Fluoro | Röntgen RF** is a customized fluoroscopy system solution. We help you creating your system. Geometries from our partners with our packages together with the development support.



## Training

- In depth product training or workshops for interfacing, basic radiology know-how or application imaging specialties; ask our engineers or book them for an in-house training.
- For all kind of technical support, service issues or spare part needs our technical support centre can be contacted.



## Online Services

- Please make use of our online services providing detailed product information, documentation, templates or other information. You can download specifications or even our general terms and conditions.
- A protected customer login zone for a close customer relationship and the need for a fast exchange of data, programs or other files.

